



BETA DIAGNOSTICI

ORANGE READER



Specifiche tecniche

Misurazioni: Qualitative
Tempo di Lettura: 10 secondi
Campioni: Urine / Feci / Siero / Plasma / Sangue Intero
Memoria Interna: 8 GB
Interfaccia di Connessione: USB (risultati esportabili in Excel o PDF)
Pistola barcode
Laptop
Software: Italiano



viale Europa, 22 is. 57 - 98123 MESSINA
Tel. +39 090 696001 Fax +39 090 695563

www.betadiagnostici.it
info@betadiagnostici.it



Qualitative Rapid Tests with Reader

Cat. No.	Product Description	Specimen	Format
IHP-R602	H. pylori Antigen Rapid Test Cassette	Feces	Cassette
IHC-R402	HCV Rapid Test Cassette	WB / S/ P	Cassette
IMO-R402	MONO Rapid Test Cassette	WB / S/ P	Cassette
IIN-R502	Influenza A+B Rapid Test Cassette	Swab / Nasal Aspirate	Cassette
IFLA-R502	Influenza A Rapid Test Cassette	Swab / Nasal Aspirate	Cassette
IST-R502	Strep A Rapid Test Cassette (Control Line in Red)	Throat Swab	Cassette
TFO-R602	FOB Rapid Test Cassette	Feces	Cassette
TPS-R402	PSA Rapid Test Cassette	WB / S/ P	Cassette
FFF-R502	Fetal Fibronectin (fFN) Rapid Test Cassette	Vaginal Secretion	Cassette
FIG-R502	Insulin-like Growth Factor-binding Protein 1 (IGFBP-1) Rapid Test Cassette	Vaginal Secretion	Cassette

Qualitative Drug of Abuse Rapid Tests with Reader

Cat. No.	Product Description	Specimen	Cut-Off
DAC-R102	Acetaminophen (ACE) Rapid Test Cassette	Urine	5000 ng/mL
DACL-R102	ACL Rapid Test Cassette	Urine	300 ng/mL
DAM-R102	Amphetamine (AMP) Rapid Test Cassette	Urine	1000 ng/mL
DBA-R102	Barbiturate (BAR) Rapid Test Cassette	Urine	300 ng/mL
DBU-R102	Buprenorphine (BUP) Rapid Test Cassette	Urine	10 ng/mL
DBZ-R102	Benzodiazepines (BZO) Rapid Test Cassette	Urine	300 ng/mL
DCL-R102	Clonazepam (CLO) Rapid Test Cassette	Urine	400 ng/mL
DCO-R102	COCAINE (COC) Rapid Test Cassette	Urine	300 ng/mL
DCT-R102	Cotinine (COT) Rapid Test Cassette	Urine	200 ng/mL
DDI-R102	Diazepam (DIA) Rapid Test Cassette	Urine	300 ng/mL
DED-R102	Ethylendiamine-dimethylphosphinic acid (EDDP) Rapid Test Cassette	Urine	100 ng/mL
DET-R102	Ethyl Glucuronide (ETG) Rapid Test Cassette	Urine	500 ng/mL
DFY-R102	Fentanyl (FYL) Rapid Test Cassette	Urine	20 ng/mL
DKE-R102	Ketamine (KET) Rapid Test Cassette	Urine	1000 ng/mL
DMA-R102	6-Monoacetylmorphine (6-MAM) Rapid Test Cassette	Urine	10 ng/mL
DMD-R102	Ecstasy (MDMA) Rapid Test Cassette	Urine	500 ng/mL
DMDA-R102	Tenamfetamine (MDA) Rapid Test Cassette	Urine	500 ng/mL
DMDP-R102	Methylenedioxypropylvalerone (MDPV) Rapid Test Cassette	Urine	1000 ng/mL
DME-R102	Methamphetamine (MET) Rapid Test Cassette	Urine	1000 ng/mL
DMO-R102	Morphine (MOP) Rapid Test Cassette	Urine	300 ng/mL
DMQ-R102	Methaqualone (MQL) Rapid Test Cassette	Urine	300 ng/mL
DMT-R102	Methadone (MTD) Rapid Test Cassette	Urine	300 ng/mL
DOP-R102	Opiates (OPI) Rapid Test Cassette	Urine	2000 ng/mL
DOX-R102	Oxycodone (OXY) Rapid Test Cassette	Urine	100 ng/mL
DPC-R102	Phencyclidine (PCP) Rapid Test Cassette	Urine	25 ng/mL
DPP-R102	Propoxyphene (PPX) Rapid Test Cassette	Urine	300 ng/mL
DSM-R102	Synthetic Marijuana (K2) Rapid Test Cassette	Urine	50 ng/mL
DTC-R102	Tricyclic Antidepressants (TCA) Rapid Test Cassette	Urine	1000 ng/mL
DAP-R102	α-Pyrrolidinovalerophenone (α-PVP) Rapid Test Cassette	Urine	1000 ng/mL
DMP-R102	Methylphenidate (MPD) Rapid Test Cassette	Urine	300 ng/mL
DAL-R102	Alcohol (ALC) Rapid Test Cassette	Urine	0.02%
DTH-R102	Marijuana (THC) Rapid Test Cassette	Urine	50 ng/mL
DTM-R102	Tramadol (TML) Rapid Test Cassette	Urine	100 ng/mL
DLS-R102	Lysergic Acid Diethylamide (LSD) Rapid Test Cassette	Urine	20 ng/mL
DZO-R102	Zolpidem (ZOL) Rapid Test Cassette	Urine	50 ng/mL
DOA-R134	Multi-Drug 3 Drugs Rapid Test Panel (MOP/AMP/THC)	Urine	N/A
DOA-R134	Multi-Drug 3 Drugs Rapid Test Panel (AMP/COC/THC)	Urine	N/A
DOA-R144	Multi-Drug 4 Drugs Rapid Test Panel (MOP/COC/THC/BZO)	Urine	N/A
DOA-R144	Multi-Drug 4 Drugs Rapid Test Panel (OPI/COC/MET/THC)	Urine	N/A
DOA-R154	Multi-Drug 5 Drugs Rapid Test Panel (MOP/COC/THC/MET/BZO)	Urine	N/A

Qualitative Drug of Abuse Rapid Tests with Reader

Cat. No.	Product Description	Specimen	Cut-Off
DOA-R154	Multi-Drug 5 Drugs Rapid Test Panel (OPI/MDMA/COC/AMP/MET)	Urine	N/A
DOA-R164	Multi-Drug 6 Drugs Rapid Test Panel (OPI/MDMA/COC/THC/AMP/MET)	Urine	N/A
DOA-R164	Multi-Drug 6 Drugs Rapid Test Panel (OXY/COC/MET/PPX/BZO/TCA)	Urine	N/A
DOA-R174	Multi-Drug 7 Drugs Rapid Test Panel (OPI/COC/THC/AMP/MET/MTD/BZO)	Urine	N/A
DOA-R174	Multi-Drug 7 Drugs Rapid Test Panel (MOP/MDMA/COC/THC/AMP/MET/PCP)	Urine	N/A
DOA-R184	Multi-Drug 8 Drugs Rapid Test Panel (MOP/COC/THC/AMP/MTD/BAR/BZO/BUP)	Urine	N/A
DOA-R184	Multi-Drug 8 Drugs Rapid Test Panel (OPI/OXY/THC/AMP/MTD/BZO/FYL/BUP)	Urine	N/A
DOA-R194	Multi-Drug 9 Drugs Rapid Test Panel (OPI/MDMA/COC/THC/AMP/MET/MTD/BZO/BUP)	Urine	N/A
DOA-R194	Multi-Drug 9 Drugs Rapid Test Panel (COC/AMP/MET/BAR/BZO/COT/FYL/TCA/BUP)	Urine	N/A
DOA-R1104	Multi-Drug 10 Drugs Rapid Test Panel (MOP/MDMA/COC/THC/AMP/MET/MTD/BAR/BZO/BUP)	Urine	N/A
DOA-R1104	Multi-Drug 10 Drugs Rapid Test Panel (OPI/OXY/COC/THC/AMP/MET/TML/BZO/FYL/TCA)	Urine	N/A
DOA-R1124	Multi-Drug 12 Drugs Rapid Test Panel (MOP/MDMA/COC/THC/AMP/MET/MTD/PCP/BAR/BZO/TCA/BUP)	Urine	N/A
DOA-R1124	Multi-Drug 12 Drugs Rapid Test Panel (OPI/COC/THC/AMP/MET/K2/MTD/PCP/BAR/BZO/TCA/ACE)	Urine	N/A
DOA-R135	Multi-Drug 3 Drugs Rapid Test Cassette (AMP/COC/THC)	Urine	N/A
DOA-R135	Multi-Drug 3 Drugs Rapid Test Cassette (MOP/AMP/MET)	Urine	N/A
DOA-R155	Multi-Drug 5 Drugs Rapid Test Cassette (AMP/COC/MET/MOP/THC)	Urine	N/A
DOA-R165	Multi-Drug 6 Drugs Rapid Test Cassette (<AMP/COC/THC>+<MET/MOP>+<BZO>)	Urine	N/A
DOA-R165	Multi-Drug 6 Drugs Rapid Test Cassette (<COC/MET/BZO>+<MDMA/MTD>+<THC>)	Urine	N/A
DOA-R175	Multi-Drug 7 Drugs Rapid Test Cassette (<AMP/COC/THC>+<MET/MOP>+<MDMA/MTD>)	Urine	N/A
DOA-R175	Multi-Drug 7 Drugs Rapid Test Cassette (<MET/COC/THC>+<AMP/OPI>+<TCA/BAR>)	Urine	N/A
DOA-R175	Multi-Drug 7 Drugs Rapid Test Cassette (<COC/MET/BZO>+<AMP/MOP>+<MDMA/MTD>)	Urine	N/A
DOA-R1125	Multi-Drug 12 Drugs Rapid Test Cassette (<MET/COC/THC>+<MDMA/MTD/MOP>+<PCP/BZO/AMP>+<TCA/BAR>+<BUP>)	Urine	N/A
DOA-R1125	Multi-Drug 12 Drugs Rapid Test Cassette (<COC/MET/BZO>+<THC/AMP/OPI>+<MDMA/MTD/MOP>+<TCA/BAR>+<ETG>)	Urine	N/A